

**Item Name :** Bernoulli's Tube

**Item Code :** 101A.2008.01

**Item Description :**

Description

To demonstrate the Bernoulli's effect produced by water flow. Comprises a Bernoulli's tube, with central constriction a constant bore gradient tube, 6 lengths glass tubing, 2 swan necked outlets tubes and 8 precut lengths of rubber connecting tubing. The 2 flow tubes, each 500x13mm, length 25mm x 150mm spacing extended to form 500mm manometers with the lengths of glass tubing.

**Item Image :**



Cat. No.

101A.2008.01