

Item Name : Burette

Item Code : 804.2008.01

Item Description :

This Burette has a transparent acrylic body and a leakproof stopcock made fitted with a self lubricating PTFE plug. Each Burette is individually calibrated to give maximum accuracy. Acrylic body is unaffected by dilute mineral acids & bases except Hydrofloric Acid & Ammonium Hydroxide. The pH range of 1-14 is ideal and can be used up to 5% wt/wt of acids or 30 % bases. Material : Polymethyl Methacrylate /TPX

Item Image :



Cat. No.	Particulars	Packing
804.2008.01	25 ml	10 Pcs
804.2008.02	50ml	10 Pcs
804.2008.03	100 ml	10 Pcs