

Item Name : Spare Burette length without stopcock

Item Code : 125.2008.01

Item Description :

class B complies with BS 846

Item Image :



Cat. No.	Cap. ml.	Class	Sub. Div. ml.	Tolerance (\pm ml)
125.2008.01	10	B	0.05	0.05
125.2008.02	25	B	0.10	0.10
125.2008.03	50	B	0.10	0.10
125.2008.04	100	B	0.20	0.20
125.2008.05	10	A	0.05	0.02
125.2008.06	25	A	0.10	0.05
125.2008.07	50	A	0.10	0.05
125.2008.08	100	A	0.20	0.10



44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>