

**Item Name :** Burette Automatic Zero

**Item Code :** 130.2008.01

**Item Description :**

Mounted on Reservoir with screw type PTFE Rotaflow Needle Valve Stopcock Schellback and Rubber bellow, Accuracy as per class 'A' o

**Item Image :**



| Cat. No.    | Cap. ml. | Sub. Div. ml. | Tolerance ( $\pm$ ml) |
|-------------|----------|---------------|-----------------------|
| 130.2008.01 | 10       | 0.05          | 0.02                  |
| 130.2008.02 | 25       | 0.10          | 0.05                  |
| 130.2008.03 | 50       | 0.10          | 0.05                  |
| 130.2008.04 | 100      | 0.20          | 0.10                  |