

Item Name : Burette Automatic Zero

Item Code : 130.2008.01

Item Description :

Mounted on Reservoir with screw type PTFE Rotaflow Needle Valve Stopcock Schellback and Rubber bellow, Accuracy as per class 'A'

Item Image :



Cat. No.	Cap. ml.	Sub. Div. ml.	Tolerance (\pm ml)
130.2008.01	10	0.05	0.02
130.2008.02	25	0.10	0.05
130.2008.03	50	0.10	0.05
130.2008.04	100	0.20	0.10