

Item Name : Burette Automatic Zeroo

Item Code : 131.2008.01

Item Description :

Mounted on Reservoir with PTFE Rotaflow Stopcock, Schellback Rubber bellow, Accuracy as per class B Dissected in 3 parts

Note : Burette Automatic Zero can also be supplied without schellback in Glass Stopcock and in P.T.F.E. Stopcock on request.

Item Image :



Cat. No.	Cap. ml.	Sub. Div. ml.	Tolerance (\pm ml)
131.2008.01	10	0.05	0.05
131.2008.02	25	0.10	0.10
131.2008.03	50	0.10	0.10
131.2008.04	100	0.20	0.20