

Item Name : Graduated Pipettes Serological

Item Code : 135.2008.01

Item Description :

Class B ISO 835, BS 700

Note :Pipette Mohr Type can also be supplied on request.

Item Image :



Cat. No.	Cap. ml.	Tolerance (\pm ml)	Colour code
135.2008.01	1.0	0.010	Yellow
135.2008.02	2.0	0.020	Black
135.2008.03	5.0	0.050	Red
135.2008.04	10.0	0.100	Orange
135.2008.06	25.0	0.200	White
135.2008.07	50.0	-	Blue
135.2008.08	100.0	-	White

A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>