

Item Name : Graduated Pipette Serological

Item Code : 136.2008.01

Item Description :

Class A, ISO835, BS 700

Note : Pipette Mohr Type can also be supplied on request.

For Batch Certificate Put B in front of Cat No.

For Individual Certificate Put C in front of Cat No.

Item Image :



Cat. No.	Cap. ml.	Tolerance (\pm ml)	Colour code
136.2008.01	1.0	0.006	Yellow
136.2008.02	2.0	0.010	Black
136.2008.03	5.0	0.030	Red
136.2008.04	10.0	0.050	Orange
136.2008.06	25.0	0.100	White
136.2008.07	50.0	-	Blue
136.2008.08	100.0	-	White

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>