

Item Name : Volumetric Flask

Item Code : 888.2008.01

Item Description :

Glacier Volumetric Flasks, blown in Polypropylene, are strong, unbreakable and have good contact clarity. These autoclavable flasks are provided with screw caps to prevent leakage. Material :PP

Item Image :



Cat. No.	Particulars	Packing
888.2008.01	25 ml	6 Pcs
888.2008.02	50 ml	6 Pcs
888.2008.03	100 ml	6 Pcs
888.2008.04	250 ml	3 Pcs
888.2008.05	500 ml	3 Pcs
888.2008.06	1000 ml	2 Pcs

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>