

Item Name : Volumetric Flask

Item Code : 888.2008.01

Item Description :

Glacier Volumetric Flasks, blown in Polypropylene, are strong, unbreakable and have good contact clarity. These autoclavable flasks are provided with screw caps to prevent leakage. Material :PP

Item Image :



| Cat. No. | Particulars | Packing |
|-------------|-------------|---------|
| 888.2008.01 | 25 ml | 6 Pcs |
| 888.2008.02 | 50 ml | 6 Pcs |
| 888.2008.03 | 100 ml | 6 Pcs |
| 888.2008.04 | 250 ml | 3 Pcs |
| 888.2008.05 | 500 ml | 3 Pcs |
| 888.2008.06 | 1000 ml | 2 Pcs |

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>