

**Item Name :** Draining Rack

**Item Code :** 933.2008.01

**Item Description :**

Draining Rack provides an excellent option for post wash draining of tubes & bottles in the laboratory. The back plate has 20 holes into which pegs are tightened with the help of screw-nuts. Bottom of the rack is provided with a tray which collects all the waste water drained down which can further be collected in a container via drain tube. This rack does not rot, stain or rust. Material : Acrylonitrile Butadiene Styrene/Polycarbonate

**Item Image :**



Cat. No.	Particulars	Packing
933.2008.01	Draining Rack	12 Pcs.