

**Item Name :** Instrument Tray

**Item Code :** 932.2008.01

**Item Description :**

This autoclavable trays have cover which fits into groove of trays to minimize spillage of material. Smaller size can be used to sterilize small instruments whereas the larger size allows pipettes to be fully immersed in disinfectant prior to sterilizing. Material : Polypropylene

**Item Image :**



Cat. No.	Particulars	Packing
932.2008.01	220 x 150 x 70mm	2 Pcs
932.2008.02	450 x 150 x 70 mm	2 Pcs