

Item Name : Volumetric Flask Class A, One Graduation Mark

Item Code : 141.2008.01

Item Description :

Stopper made of polythelene ISO 1042,DIN 12664

Item Image :



Cat. No.	Cap. ml.	Tolerance (\pm ml)	N/S
141.2008.01	5	0.025	10/19
141.2008.02	10	0.025	10/19
141.2008.02A	20	0.040	10/19
141.2008.03	25	0.040	10/19
141.2008.04	50	0.060	12/21
141.2008.05	100	0.100	14/23
141.2008.06	200	0.150	14/23
141.2008.07	250	0.150	14/23
141.2008.08	500	0.250	19/26
141.2008.09	1000	0.400	19/26

141.2008.10	2000	0.600	29/32
141.2008.11	5000	1.200	34/35

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>