

**Item Name :** Laboratory Tray

**Item Code :** 929.2008.01

**Item Description :**

Laboratory Tray moulded in polypropylene is a multipurpose tray that can be used for sterilizing, drying glassware, porcelain ware & other laboratory instruments. These steam autoclavable trays have tapered walls which prove good for nesting. Material : Polypropylene

**Item Image :**



Cat. No.	Particulars	Packing
929.2008.01	450 x 350 x 75 mm	10 Pcs
929.2008.02	375 x 300 x 75 mm	10 Pcs