

Item Name : Plantation Pots

Item Code : 928.2008.01

Item Description :

Moulded in polypropylene, are available in four different sizes. These pots have a downward extended rim which provide convenient lifting even if the pot is filled. the bottom of the pot is provided with holes for drain purpose. Material : Polypropylene

Item Image :



| Cat. No. | Particulars | Packing |
|-------------|-------------|---------|
| 928.2008.01 | 100 mm | 20 Pcs |
| 928.2008.02 | 150 mm | 20 Pcs |
| 928.2008.03 | 200 mm | 20 Pcs |
| 928.2008.04 | 250 mm | 20 Pcs |