

Item Name : Calorimeter With Wooden Jacket (Stirrer + Lid Superior)

Item Code : 3014.2008.01

Item Description :

Wooden Joule Calorimeter, The size of this calorimeter is 3"X2"

Following are the features of wooden calorimeter

- Electrically heated calorimeter, such as that used in the Griffiths method.
- It consists of a copper vessel with outer vessel as a superior quality teak wood box of size 9.5 X 9.5 X 11.5 cm fitted with Bakelite lid with hole for the thermometer and stirrer.



Item Image :

Cat No.

3014.2008.01