

Item Name : Volumetric Flask Class A, One Graduation Markk

Item Code : 141.2009.01

Item Description :

Stopper made of polythelene ISO 1042, DIN 12664 with Batch Certificate

Item Image :



Cat. No.	Cap. ml.	Tolerance (\pm ml)	N/S
141.2009.01	5	0.025	10/19
141.2009.02	10	0.025	10/19
141.2009.02A	20	0.040	10/19
141.2009.03	25	0.040	10/19
141.2009.04	50	0.060	12/21
141.2009.05	100	0.100	14/23
141.2009.06	200	0.150	14/23
141.2009.07	250	0.150	14/23
141.2009.08	500	0.250	19/26
141.2009.09	1000	0.400	19/26

141.2009.10	2000	0.600	29/32
141.2009.11	5000	1.200	34/35

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>