


**Item Name :** Volumetric Flask Class A, One Graduation Marrk

**Item Code :** 141.2010.01

**Item Description :**

Stopper made of polythelene ISO 1042, DIN 12664 with Individual Works Certificate  
 Serially Numbered

**Item Image :**



Cat. No.	Cap. ml.	Tolerance ( $\pm$ ml)	N/S
141.2010.01	5	0.025	10/19
141.2010.02	10	0.025	10/19
141.2010.02A	20	0.040	10/19
141.2010.03	25	0.040	10/19
141.2010.04	50	0.060	12/21
141.2010.05	100	0.100	14/23
141.2010.06	200	0.150	14/23
141.2010.07	250	0.150	14/23
141.2010.08	500	0.250	19/26
141.2010.09	1000	0.400	19/26

141.2010.10	2000	0.600	29/32
141.2010.11	5000	1.200	34/35

# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

[info@glacierindia.com](mailto:info@glacierindia.com)

<https://www.glacierindia.com>