

Item Name : Volumetric Flask Amber, one Graduation mark

Item Code : 142.2008.01

Item Description :

Stopper made of Polythelene ISO 1042 DIN 12664. Accuracy as per class B

Item Image :



Cat. No.	Cap. ml.	Tolerance (±ml)	N/S
142.2008.01	5	0.05	10/19
142.2008.02	10	0.05	10/19
142.2008.03	20	0.08	10/19
142.2008.04	25	0.08	10/19
142.2008.05	50	0.12	12/21
142.2008.06	100	0.20	14/23
142.2008.07	200	0.30	14/23
142.2008.08	250	0.30	14/23
142.2008.09	500	0.50	19/26
142.2008.10	1000	0.80	19/26

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>