


**Item Name :** Universal Reagent Reservoir

**Item Code :** 908.2008.01

**Item Description :**

This Universal Reagent Reservoir is a graduated single reagent reservoir with maximum volume of 50 ml, reservoir with 8 separate numbered 5ml. reagent channels & a reservoir with 12 separate numbered 5 ml. reagent channels. Made of polypropylene these reagent reservoirs are repeatedly autoclavable and can be used with all multi-channel pipettes. Material : Polypropylene

**Item Image :**



Cat. No.	Particulars	Packing
908.2008.01	50 ml	6 Pcs