

Item Name : Volumetric Flask Amber Class A

Item Code : 143.2009.01

Item Description :

One Graduation mark & stopper made of Polythelene ISO 1042 DIN 12664. with Batch Certificate

Item Image :



Cat. No.	Cap. ml.	Tolerance (±ml)	N/S
143.2009.01	5	0.025	10/19
143.2009.02	10	0.025	10/19
143.2009.03	20	0.040	10/19
143.2009.04	25	0.040	10/19
143.2009.05	50	0.060	12/21
143.2009.06	100	0.100	14/23
143.2009.07	200	0.150	14/23
143.2009.08	250	0.150	14/23
143.2009.09	500	0.250	19/26
143.2009.10	1000	0.400	19/26

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>