

Item Name : Potometerr

Item Code : 146.2008.01

Item Description :

Ideal for laboratory & field use, it allows the students to measure transpiration with minimal risk of breakages. The apparatus has only one glass component and is mounted on a strong backing plate. It consists of a capillary tube with mm scale connected via a 3 way tap to a syringe. The cut shoot is connected to the third outlet. The apparatus can also be used to measure humidity by connecting a porous pot to the third outlet. Capillary tube 120 x 1mm length x Bore

Item Image :



Cat. No.

146.2008.01