

**Item Name :** Test Tube stand 3 Tier

**Item Code :** 897.2008.01

**Item Description :**

These ingenious racks with a 3-tier design provide clear view of the tube contents & convenient insertion & removal of tubes. The PC rack has transparent upper & middle plates where as in PP racks, every part is coloured & opaque. These autoclavable Racks can withstand sub- freezing temperatures. Material : Polycarbonate

**Item Image :**



Cat. No.	Particulars	Packing
897.2008.01	13 mm x 31 Tubes	2 Pcs
897.2008.02	16 mm x 31 Tubes	2 Pcs
897.2008.03	20 mm x 20 Tubes	2 Pcs
897.2008.04	25 mm x 12 Tubes	2 Pcs
897.2008.05	13 mm x 62 Tubes	2 Pcs
897.2008.06	16 mm x 62 Tubes	2 Pcs
897.2008.07	20 mm x 40 Tubes	2 Pcs
897.2008.08	25 mm x 24 Tubes	2 Pcs
897.2008.09	25 mm x 18 Tubes	2 Pcs

897.2008.10      25 mm x 36 Tubes      2 Pcs

897.2008.11      32 mm x 12 Tubes      2 Pcs

# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

[info@glacierindia.com](mailto:info@glacierindia.com)

<https://www.glacierindia.com>