


Item Name : Test Tube Stand Wire Pattern

Item Code : 895.2008.01

Item Description :

An economical substitute for wire racks, here are these submersible, autoclavable racks which are stackable when empty. These are available in different colours and are made of a special blend of Polypropylene that makes it sink in a water bath and maintain their stability even in agitated water. These space saving racks can easily be disassembled for cleaning purpose. The base is in the form of a grid with square openings and has fold up sides and end bars that lock in to rack tops. Material : Polypropylene

Item Image :



Cat. No.	Particulars	Packing
895.2008.01	13 mm x 90 Tubes	6 Pcs
895.2008.02	16 mm x 60 Tubes	6 Pcs
895.2008.03	20 mm x 40 Tubes	6 Pcs
895.2008.04	25 mm x 24 Tubes	6 Pcs
895.2008.05	30 mm x 21 Tubes	6 Pcs

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>