

Item Name : Test Tube stand 3 Tier PP

Item Code : 894.2008.01

Item Description :

These ingenious racks with a 3-tier design provide clear view of the tube contents & convenient insertion & removal of tubes. The PC rack has transparent upper & middle plates where as in PP racks, every part is colored & opaque. These autoclavable Racks can withstand sub- freezing temperatures. Material : Polypropylene



Item Image :

Cat. No.	Particulars	Packing
894.2008.01	13 mm x 31 Tubes	2 Pcs
894.2008.02	16 mm x 31 Tubes	2 Pcs
894.2008.03	20 mm x 20 Tubes	2 Pcs
894.2008.04	25 mm x 12 Tubes	2 Pcs
894.2008.05	25 mm x 18 Tubes	2 Pcs
894.2008.06	32 mm x 12 Tubes	2 Pcs
894.2008.07	25 mm x 36 Tubes	2 Pcs
894.2008.08	13 mm x 18 Tubes	2 Pcs
894.2008.09	15 mm x 18 Tubes	2 Pcs

894.2008.10	18 mm x 12 Tubes	2 Pcs
894.2008.11	13 mm x 62 Tubes	2 Pcs
894.2008.12	16 mm x 62 Tubes	2 Pcs
894.2008.13	20 mm x 40 Tubes	2 Pcs
894.2008.14	13 mm x 48 Tubes	2 Pcs
894.2008.15	13 mm x 72 Tubes	2 Pcs

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>