

Item Name : Measuring Cylinders Plastic Bases

Item Code : 151.2008.01

Item Description :

Protection collar Class B

Item Image :



Cat. No.	Cap. ml.	Sub. Div. ml.	Tolerance (\pm ml)
151.2008.01	5	0.1	0.1
151.2008.02	10	0.2	0.2
151.2008.03	25	0.5	0.5
151.2008.04	50	1.0	1.0
151.2008.05	100	1.0	1.0
151.2008.06	250	2.0	2.0
151.2008.07	500	5.0	5.0
151.2008.08	1000	10.0	10.0
151.2008.09	2000	20.0	20.0

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>