

**Item Name :** Test Tube Stand Round

**Item Code :** 893.2008.01

**Item Description :**

A handy, little, space saver, this circular twelve place Polypropylene Test Tube Stand has four 25 mm. dia holes & eight 19 mm dia holes in its top plate. The base has matching hemispherical wells & 16 vertical pins for drying tubes surrounded by a trough to catch drippings. This autoclavable rack can be easily disassembled for cleaning purpose. Material : Polypropylene

**Item Image :**



Cat. No.	Particulars	Packing
893.2008.01	19 mm & 25 mm Ø Tubes	2 Pcs.