

Item Name : Test Tube Stand

Item Code : 892.2008.01

Item Description :

This autoclavable stand is perhaps the most commonly used Test Tube Stand in various laboratories. For better visibility, this 6 place rack is end supported by 2 columns. The top plate has 2 holes for holding 25 mm dia test tubes, four holes to hold 16 mm dia tubes where as the base has corresponding hemispherical wells. It also has six vertical pins for drying tubes.

Material : Polypropylene

Item Image :



Cat. No.	Particulars	Packing
892.2008.01	16 mm & 25 mm Ø Tubes	12 Pcs.
892.2008.02	25 mm Ø Tubes	12 Pcs.