

**Item Name :** Measuring Cylinders Spout & Round Base

**Item Code :** 152.2008.01

**Item Description :**

Class ADIN 12680, ISO 4788 BS 604

**Item Image :**



Cat. No.	Cap. ml.	Sub. Div. ml.	Tolerance ( $\pm$ ml)
152.2008.01	5	0.1	0.05
152.2008.02	10	0.2	0.10
152.2008.03	25	0.5	0.25
152.2008.04	50	0.1	0.50
152.2008.05	100	0.1	0.50
152.2008.06	250	2.0	1.00
152.2008.07	500	5.0	2.50
152.2008.08	1000	10.0	5.00
152.2008.09	2000	20.0	10.0

# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>