

Item Name : Micro Tip Box

Item Code : 886.2008.01

Item Description :

These unique boxes are moulded in Polypropylene and are therefore autoclavable. The Box itself is clear providing an excellent view of the tips inside while a table like structure is provided inside to hold the tips. The tip holding spaces are uniformly placed to facilitate the use of multi channel pipettes. The cover is joined to the base through ball catcher hinges and provides excellent protection to the tips. Material : Polypropylene

Item Image :



Cat. No.	Particulars	Packing
886.2008.01	For 96 Micro Tips of 2-10 ul	6 Pcs
886.2008.02	For 100 Micro Tips of 200-1000 ul	6 Pcs