

**Item Name :** Micro Test Plates

**Item Code :** 885.2008.01

**Item Description :**

ELISA/ASSAY plates are made of optically clear, non- toxic Polystyrene. These plates with 96 U Bottom or flat bottom wells bear alpha- numeric index. Uniform placement of all the wells facilitate use of multi channel pipettes. These plates packed individually are available non- sterile as well as radiation sterilized. Material : Polystyrene

**Item Image :**



Cat. No.	Particulars	Packing
885.2008.01	Flat Bottom	72 Pcs.
885.2008.02	U Bottom	72 Pcs.