


Item Name : Cryo Box PP

Item Code : 884.2008.01

Item Description :

Moulded in polypropylene, the top of the box is attached to the bottom through well designed built in hinges. The bottom has built- in spaces to hold vials up to 2 ml. capacity. Available in two different sizes, these autoclavable boxes bear alpha numeric index for easy identification of different vials placed inside. Material : Polypropylene

Item Image :



Cat. No.	Particulars	Packing
884.2008.01	81 places for 1 ml. or 1.8 ml Cryo Vials	4 Pcs
884.2008.02	100 places for 1 ml.or 1.8 ml Cryo Vials	4 Pcs