

Item Name : Measuring Cylinders Spout & Round 'Base'

Item Code : 152.2010.01

Item Description :

Class ADIN 12680, ISO 4788, BS 604, with Individual works certificate

Item Image :



| Cat. No. | Cap. ml. | Sub. Div. ml. | Tolerance (\pm ml) |
|-------------|----------|---------------|-----------------------|
| 152.2010.01 | 5 | 0.1 | 0.05 |
| 152.2010.02 | 10 | 0.2 | 0.05 |
| 152.2010.03 | 25 | 0.5 | 0.25 |
| 152.2010.04 | 50 | 1.0 | 0.50 |
| 152.2010.05 | 100 | 1.0 | 0.50 |
| 152.2010.06 | 250 | 2.0 | 1.00 |
| 152.2010.07 | 500 | 5.0 | 2.50 |
| 152.2010.08 | 1000 | 10.0 | 5.00 |
| 152.2010.09 | 2000 | 20.0 | 10.0 |

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>