

Item Name : Measuring Cylinders Spout & Round 'Base'

Item Code : 152.2010.01

Item Description :

Class ADIN 12680, ISO 4788, BS 604, with Individual works certificate

Item Image :



Cat. No.	Cap. ml.	Sub. Div. ml.	Tolerance (\pm ml)
152.2010.01	5	0.1	0.05
152.2010.02	10	0.2	0.05
152.2010.03	25	0.5	0.25
152.2010.04	50	1.0	0.50
152.2010.05	100	1.0	0.50
152.2010.06	250	2.0	1.00
152.2010.07	500	5.0	2.50
152.2010.08	1000	10.0	5.00
152.2010.09	2000	20.0	10.0

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>