


**Item Name :** Cryo Box PC

**Item Code :** 883.2008.01

**Item Description :**

Cryo Boxes are moulded in Polycarbonate making it extra strong and autoclavable. These can accommodate Cryo Vials of up to 4.5 ml. capacity and are designed to be used at temperatures ranging from -190°C to +121°C. The transparent cover of the box gives a good view of the contents inside and is numbered for better inventory control. Material : Polycarbonate

**Item Image :**



Cat. No.	Particulars	Packing
883.2008.01	50 places for 1 ml. or 1.8 ml. Cryo Vials	8 Pcs.
883.2008.02	81 places for 1 ml. or 1.8 ml. Cryo Vials	4 Pcs.
883.2008.03	81 places for 4.5 ml. Cryo Vials	4 Pcs.
883.2008.04	100 places for 1 ml. or 1.8 ml. Cryo Vials	4 Pcs.