

Item Name : Cryo Rack

Item Code : 882.2008.01

Item Description :

Moulded in Polycarbonate making it tough, durable & auto clavable . This handy rack can hold 50 vials at a time. The unique feature of this rack is that vials can be unscrewed without actually picking them up as the vial will lock in every well and hence will not twist. This rack is provided with four anti skid rubbers & each well is identified with an alpha- numeric index.

Material : Polycarbonate

Item Image :



| Cat. No. | Particulars | Packing |
|-------------|------------------------------|---------|
| 882.2008.01 | 50 places for 1,1.8 & 4.5 ml | 2 Pcs. |