

**Item Name :** Cryo Vial

**Item Code :** 880.2008.01

**Item Description :**

Cryo Vials are normally used for storage of Biological material, human & animal cells at temperature as low as  $-190^{\circ}\text{C}$ . These sterilized tubes are provided with white, prominent graduations and a long marking area for sample identification. The cap of these vials is provided with a seal ring inside to make it completely leakproof. The longitudinal grooves in the base give a non-twisting feature to these vials when placed in a cryo rack. Material : Polypropylene

**Item Image :**



Cat. No.	Particulars	Packing
880.2008.01	1.0 ml	1000 Pcs
880.2008.02	1.8 ml	1000 Pcs
880.2008.03	4.5 ml	500 Pcs