

Item Name : Cryo Vial

Item Code : 880.2008.01

Item Description :

Cryo Vials are normally used for storage of Biological material, human & animal cells at temperature as low as -190°C . These sterilized tubes are provided with white, prominent graduations and a long marking area for sample identification. The cap of these vials is provided with a seal ring inside to make it completely leakproof. The longitudinal grooves in the base give a non-twisting feature to these vials when placed in a cryo rack. Material : Polypropylene

Item Image :



Cat. No.	Particulars	Packing
880.2008.01	1.0 ml	1000 Pcs
880.2008.02	1.8 ml	1000 Pcs
880.2008.03	4.5 ml	500 Pcs