


**Item Name :** Rack for Scientillation Vial

**Item Code :** 878.2008.01

**Item Description :**

Moulded in Polypropylene, are autoclavable & corrosion resistant. There is an individual rack for both 8 ml. & 20 ml. vials. These rugged & stable racks can hold scintillation vials for convenient storage, handling & sample preparations. Material : Polypropylene

**Item Image :**



Cat. No.	Particulars	Packing
878.2008.01	Rack for 90 vials of 8ml. capacity	2 Pcs
878.2008.02	Rack for 50 vials of 20ml. capacity	2 Pcs