

**Item Name :** Storage Vial

**Item Code :** 873.2008.01

**Item Description :**

These Polypropylene / Polystyrene tubes are generally used in RIA, Coagulation & Bacteriology. The Polypropylene tubes are clear and can withstand 3000 RCF during Centrifugation while the transparent polystyrene tubes can withstand 1400 RCF. The PP tubes can accept common acid, alkalies & solvents at room temperature while the PS tubes can accept mild bases, weak acids but not organic solvents. Material : Polyethylene

**Item Image :**



Cat. No.	Particulars	Packing
873.2008.01	1.0 ml	500 Pcs
873.2008.02	2.0 ml	500 Pcs
873.2008.03	5.0 ml	500 Pcs