

**Item Name :** Ria Vial

**Item Code :** 870.2008.01

**Item Description :**

These Polypropylene / Polystyrene tubes are generally used in RIA, Coagulation & Bacteriology. The Polypropylene tubes are clear and can withstand 3000 RCF during Centrifugation while the transparent polystyrene tubes can withstand 1400 RCF. The PP tubes can accept common acid, alkalies & solvents at room temperature while the PS tubes can accept mild bases, weak acids but not organic solvents. Material : Polypropylene / Polystyrene

**Item Image :**



Cat. No.	Particulars	Packing
870.2008.01	12 x 75 mm / PS	500 Pcs.
870.2008.02	12 x 75 mm / PP	500 Pcs.
870.2008.03	Cap for 12 mm RIA Tube	100 Pcs.