

**Item Name :** Oak Ridge Centrifuge Tubes

**Item Code :** 866.2008.01

**Item Description :**

Oak Ridge Centrifuge Tubes made of Polypropylene, have excellent contact clarity and can be used upto 50,000xG. The screw cap is fitted with silicon Oo Ring making it totally leakproof. These tubes should not be used below 4 C. Material : Polypropylene

**Item Image :**



Cat. No.	Particulars	Packing
866.2008.01	30 ml	12 Pcs
866.2008.02	50 ml	12 Pcs
866.2008.03	100 ml	12 Pcs