

**Item Name :** Centrifuge Tube Round Bottom

**Item Code :** 865.2008.01

**Item Description :**

Centrifuge tube with round bottom are moulded in Polypropylene & Polycarbonate These tubes have excellent clarity and are strong enough to be used up to 50,000x G. These tubes can be autoclaved, however, repeated autoclaving in PC tubes reduces their mechanical. Material : Polycarbonate / Polypropylene

**Item Image :**



Cat. No.	Particulars	Packing
865.2008.01	16 ml. cent. Tube PP	50 Pcs
865.2008.02	50 ml. cent. Tube PP	50 Pcs
865.2008.03	16 ml. cent. Tube PC	25 Pcs
865.2008.04	50 ml. cent. Tube PC	25 Pcs
865.2008.05	Cap for 16 ml. Cent. Tube	100 Pcs
865.2008.06	Cap for 50 ml. Cent. Tube	100 Pcs

A/Cantt-133 006. Haryana, India.

+0171 – 2698107

[info@glacierindia.com](mailto:info@glacierindia.com)

<https://www.glacierindia.com>