

**Item Name :** Centrifuge Tube Conical Bottom

**Item Code :** 864.2008.01

**Item Description :**

Centrifuge tube with conical bottom are available in two sizes i.e. 15 ml. & 50 ml. Moulded in Polypropylene, these tubes are clear, strong & safe. These autoclavable tubes are provided with printed graduations on the outer surface. The built in seal ring in the cap makes these tubes absolutely leak proof. These tubes can withstand forces up to 7500 x G & 12000 x G for 15 ml. & 50 ml. tubes. Material : Polypropylene+

**Item Image :**



Cat. No.	Particulars	Packing
864.2008.01	15 ml Non Sterile	100 Pcs
864.2008.02	50 ml Non Sterile	100 Pcs
864.2008.03	15 ml-sterile	100 Pcs
864.2008.04	50 ml-sterile	100 Pcs