

Item Name : Measuring Cylinders Round Base

Item Code : 158.2009.01

Item Description :

Interchangeable Polythelene stopper, Class A DIN 12685, ISO 4788 BS 604, with Batch Certificate

Note : Glass Stopper can be supplied on request at extra cost.

Item Image :



Cat. No.	Cap. ml.	Sub. Div. ml.	Tolerance (\pm ml)
158.2009.01	10	0.2	0.10
158.2009.02	25	0.5	0.25
158.2009.03	50	1.0	0.50
158.2009.04	100	1.0	0.50
158.2009.05	250	2.0	1.00
158.2009.06	500	5.0	2.50
158.2009.07	1000	10.0	5.00
158.2009.08	2000	10.0	10.0

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>