

Item Name : Measuring Cylinders Round Basse

Item Code : 158.2008.01A

Item Description :

Interchangeable Polythelene stopper, Class A DIN 12685, ISO 4788, BS 604, with Individual works certificate

Note : Glass Stopper can be supplied on request at extra cost.

Item Image :



Cat. No.	Cap. ml.	Sub. Div. ml.	Tolerance (\pm ml)
158.2008.01A	10	0.2	0.10
158.2008.02A	25	0.5	0.25
158.2008.03A	50	1.0	0.50
158.2008.04A	100	1.0	0.50
158.2008.05A	250	2.0	1.00
158.2008.06A	500	5.0	2.50
158.2008.07A	1000	10.0	5.00
158.2008.08A	2000	10.0	10.0

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>