

**Item Name :** Measuring Cylinders Round Base

**Item Code :** 159.2008.01

**Item Description :**

Interchangeable Polythelene stopper, Class B, DIN 12685, ISO 4788 BS 604

**Item Image :**



| Cat. No.    | Cap. ml. | Sub. Div. ml. | Tolerance ( $\pm$ ml) |
|-------------|----------|---------------|-----------------------|
| 159.2008.01 | 10       | 0.2           | 0.2                   |
| 159.2008.02 | 25       | 0.5           | 0.5                   |
| 159.2008.03 | 50       | 1.0           | 1.0                   |
| 159.2008.04 | 100      | 1.0           | 1.0                   |
| 159.2008.05 | 250      | 2.0           | 2.0                   |
| 159.2008.06 | 500      | 5.0           | 5.0                   |
| 159.2008.07 | 1000     | 10.0          | 10.0                  |
| 159.2008.08 | 2000     | 20.0          | 20.0                  |



# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>