

Item Name : Nessler Cylinder Class A, Graduated

Item Code : 162.2008.01

Item Description :

for colour comparison .

Item Image :



Cat. No.	Cap. ml.	Sub. Div. ml.	Tolerance (\pm ml)
162.2008.01	50	25 & 50	0.4
162.2008.02	100	50 & 100	0.8
162.2008.03	50	25 & 50	0.4
162.2008.04	100	50 & 100	0.8
162.2008.05	50	25 & 50	0.4
162.2008.06	100	50 & 100	0.8