


Item Name : Displacement Vessel

Item Code : 167.2008.01

Item Description :

For experiments in density, specific gravity etc. Where volumetric displacement of liquids is involved. Seamless fabrication with flared top edge, made of borosilicate glass with angled spout for convenient overflow and collection of liquids.

Item Image :



Cat. No.	Description	Capacity (ml)
167.2008.01	Graduated	250ml
167.2008.02	Graduated	500ml
167.2008.03	Without graduation	250ml
167.2008.04	Without graduation	500ml