

Item Name : Retorl Stand Mild Steel

Item Code : 1101.2008.01

Item Description :

Retorl Stand, Mild Steel the base is made of mild steel sheet, die pressed, finished with powder coating blue colour, with black rubber grip. Supplied with rod. As per DIN 12829

Item Image :



Cat. No.	Shape	Base size mm	Rod Length mm	Rod dia
1101.2008.01	Rectangular	160x100mm	500mm	12mm
1101.2008.02	Rectangular	200x125mm	600mm	12mm
1101.2008.03	Rectangular	250x160mm	750mm	12mm